

Phage Propagation Liquid Amplification (Large-Scale Lysate) Protocol

2026年2月28日 星期六

To amplify bacteriophage using *Xanthomonas* host bacteria and generate high-titer phage lysate. Used to detect phage and confirm infection. Used to expand phage stock.

Materials

- Log-phase *Xanthomonas* culture ($OD_{600} \approx 0.4-0.6$)
- Phage stock
- NYG soft agar (0.6–0.7%)
- NYG agar plates
- NYG liquid medium
- 30°C incubator

Procedure

1. Grow 50 mL *Xanthomonas* to $OD_{600} \sim 0.4$.
2. Add phage at $MOI \sim 0.01-0.1$.
 - a. MOI (Multiplicity of Infection) is calculated by dividing the total number of phage particles added by the total number of bacterial cells present.
$$\text{total phage} = \text{phage titer (PFU/mL)} \times \text{volume added (mL)}$$
$$\text{total bacteria} = \text{cell concentration (cells/mL)} \times \text{culture volume (mL)}$$
$$MOI = \text{total phages} \div \text{total bacteria}$$
3. Incubate at 30°C, 200 rpm.
4. Monitor culture:
 - Culture becomes clear/turbid drop → lysis occurring.
 - Typically 4–8 hours.
5. After complete lysis:
 - Centrifuge at 4000 x g for 10 min to remove debris.
 - Collect supernatant (contains phage).
6. Filter sterilize (0.22 μm filter).

This is crude phage lysate.

Phage Storage

- Store lysate at 4°C (short-term).
- Add chloroform (few drops) if long-term bacterial contamination risk.
- Avoid freezing unless protocol validated.

☒ Explore the Notebook 